

Sony UP25MD Specs

SONY



[Request Price Quote](#)

[Request Service Quote](#)

[Product Overview](#)



Sony UP25MD Color Printer Specs

Printing Method	Dye Sublimation Thermal Printing
Thermal Head	423 DPI
Gradations	YMC images using 8 bit/256 gradations per color, approx. 16.77 million colors per pixel
Resolution	423 DPI
Printed Picture Size	UPC-21L/UPC-24LA Size: Full 5" X 3 3/4" UPC-21S/UPC-24SA Size: Full 3 3/4" X 2 7/8"
Print Speed	NTSC: UPC-21S approx. 19 sec. UPC-21L approx. 29 sec.
Print Speed (Laminated Media)	NTSC: UPC-24SA approx. 25 sec. UPC-24LA approx. 36 sec.
Paper Capacity	S size tray: Max. 80 sheets L size tray: Max 50 sheets
Paper Capacity Laminated	S size tray: Max. 60 sheets L size tray: Max 40 sheets
Power Consumption	AC 100 to 240V, 50/60Hz
Input Current	1.7-1.0 A
Dimensions	8 3/8" X 3 7/8" X 15 5/8" (WHD)
Weight	12.6 lbs.
Video Input	RGB/YPbPr (BNC X 3) S Video (4-pin min-DIN X 1) Composite Video (BNC X 1) UPC-21L : Color Printing Pack (L size, 50 sheets x 4) UPC-21S : Color Printing Pack (S size, 80 sheets x 3) UPC-24LA : Self-Laminating Color Printing Pack (L size, 40 sheets x 4) UPC-24SA : Self-Laminating Color Printing Pack (S size, 60 sheets x 3)
Media	UL60601-1, UL60950-1 CSA C22.2 No.601.1, CSA C22.2 No.60950-1
Safety Standards	EMC: FCC: Part 15 Subpart B Digital Device Class A IC: ICES-003 Class A Digital Apparatus
Accessories	AC Power Cable Instruction Manual (PDF) Warranty Card (1) (for customers in the U.S.A. and Canada) Paper tray (L Size) Stopper (1) Printing Media (UPC-21L, 50 sheets & 1 ribbon) Cleaning Cartridge (1)