



9 GJ!- ') Ventilator



- Classical emergency transport ventilator
- Pneumatically driven pneumatically controlled
- With A/C mode and pressure trigger
- With pressure alarm
- Oxygen concentration adjustable
- Adult/pediatric parameter in different color, easy to operate
- Simple hose connecting, time saving
- Compact and portable
- Can be used at emergency department, emergency station, county hospital, clinic and for patient's inter-departmental transfer

Technical Specifications

Specification

- Modes of ventilation
Assist/Control (A/C)
Manual
Suction function
- Parameter
MV : 3~20 L/min (used in 60%)
3~14L/min (used in 100%)
Respiratory rate: 10~60bpm
I:E ratio: 1:2
FiO₂: 60%,100%
Trigger sensitivity: Pressure, -3cmH₂O
- Monitoring
- Pressure values: P_{aw}
- Alarm: P_{aw} high, P_{aw} low, O₂ supply down

Technical data

- Gas supply: O₂, 0.25~0.6MPa
- Maximum security pressure: 8kPa
- Compliance: 4mL/100Pa
- Noise: 65dB(A)
- Dimensions (H×W×D) : 377×110×142mm
- Weight: 2.4kg



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

