



ES-8000D

8 Y'jj Yf m7 UV]bYh



DENTISTRY



ES-120P

Delivery Cabinet



Cod. **ES-120P**

Standard configuration:

- 4-hole/2hole instrument hose - 2pcs
- 3-way syringe - 1pc
- Vacuum Package with Saliva ejector
- Self Contained Water Supply 1set
- Variable speed disc-type foot control
- Water cool and on/off switch with flow control
- A round handpiece pressure gauge
- Self Contained Oil less compressor (580W)
- 110V or 220V Compatible
- 10L tank capacity

Optional:

- Curing light
- Scaler

The ES-120P Delivery Cabinet offers a pullout breadboard, two instrument drawers, one medium storage drawers for materials and supplies that fully extend and quickly remove for easy cleaning. In addition, these cabinets provide a large open space for accessories with convenient access from the front.

The oversized dual-caster wheels, combined with the four-stationary handles, makes positioning of the cabinet simple and easy on all floor surfaces, even carpet. Six easy to remove and clean full extension drawers and one breadboard. High quality hardware throughout. Oversized caster wheels. Four stationary handles for easier cart positioning. A removable top plastic tray that is sectioned for instrument and accessory placement.



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

