



ES-PXD1

Analogic dental
x-ray panoramic unit



With Ceph



.....89BHGF M

ES-PXD1

Analogic dental x-ray panoramic unit



Model	ES-PXD1
Type	Oral high frequency panoramic X-ray machine
Tube voltage	60-88 KV
Tube current	Panoramic photography 4~10mA Skull photography 12mA
Photography mode	Manual / automatic setting
Rectifier mode	High frequency DC frequency 60~80 kHz
Power Supply	AC220V/110V 50/60Hz 2.5kVA
X-ray tube	Toshiba D-054SB
focus	0.5mm x 0.5mm
Total filtration	≥ 2.5mmAl
Magnification	1.30 - 1.25 times
Panoramic photography	Adult / child Adult panoramic 15 seconds Children's panoramic 12 seconds
Exposure Time	79D< 0.2-2.5 gYV6bXg f2.4-30mA.s) Temporomandibular joint 3x4 seconds
Temporomandibular joint photography	Four-segment photography
Lifting method	Electric lift
Dimensions	1880 x 881 x 2045 mm
Weight	180 kg