

ES-S300



ES-S300 is 300W analog model with basic cut, coag and bipolar modes with capability of doing most of the surgeries with lowest price

Quick reference


Power – 300 W
Analog
Non Programmable
Capability – basic surgeries
Unbelievable economy for basic functionality



1. High power efficiency of 97 %

The power efficiency of the generators are the highest in its class at 97 %. This leads to multiple benefits like increased reliability, most compact size, high portability, economy and reduced power consumption.

	Max. Power (Watts)	Rated Load (Ohms)
Body	plastic	-
Operating Frequency	480 KHz	-
Pure cut	300 W	500
Coagulation	120 W	500
Bipolar	80 W	500
Bipolar Coagulation Macro	80 W	500
Dimension (W x H x D)	12H x 33W x 36L cm	
Weight	3.0 Kg	
Input Voltage	170 to 260 VAC 50 Hz	



STANDARD ACCESSORIES



Disposable Hand Switch Pencil



Bipolar Forcep Autoclavable



Silicon Patient Plate



Chuck Handle with Electrodes Autoclavable



Monopolar Double Pedal Footswitch



ESSE 3 Via Garibaldi 30
 14022Castelnuovo D.B. (AT)
 tel +39 011 99 27 706
 fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web : www.4ci.com

