



# ESZC Line

## Full Automatic ESR/HCT Analyzers



LABORATORY



# ESZC Line

## Full Automatic ESR/HCT Analyzers



- 6.5 inch large screen, operating state of all channels in real-time display
- Track, scan and test functions of sample liquid surface
- Apply to vacuum blood tube and test tube
- Separate test function of ESR and HCT
- Buzzer alarm function after test
- Combined test function of ESR and HCT
- ESR dynamic curve display and print function
- Automatic conversion ESR equation K value function
- Test function of ESR Vm, Tm
- Stagnation time test function before erythrocyte sedimentation
- Query and print function for single and batch test results
- Test results report and data transmission function

Model	ESZC-30	ESZC-60	ESZC-100
Sample Sites (holes)	30	60	100
ESR Test Speed (Test/h)	60	120	200
HCT Test Speed (Test/h)	4300	8600	15000
Blood	1.6 ml	1.6 ml	1.6 ml



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

