



25585

Infrared multi-function forehead thermometer



LITTLE ACCESSORIES



25585

Infrared multi-function forehead thermometer





cod. **25585**

Multifunctional innovative infrared thermometer with 5 applications:
forehead temperature, room temperature, scan mode for surface temperatures (bath water or milk), real time clock, stopwatch to monitor the heart beat.


- New 3 LED zone indicator to get fever status friendly
- Instant measurement
- Large LCD easy to read
- Detailed 9 memories with time/date & site
- Fever alarm with beeper
- Selectable °C/°F mode
- Small size to carry conveniently when travelling
- Accurate: $\pm 0.2/0.3$ °C
- Self-diagnosis warning when device is abnormal
- Packed in deluxe carton box

5 APPLICATIONS



Technical Specifications:

- Measurement Range:
 - forehead 34,0°C~42,2°C
93,2°F~108,0°
 - scan mode: -22/80°C
-7,6-/76°F
- Response time: 5-30 sec
- Time between Measurement: < 5 sec
- Display Model ORAL °C/°F
- Battery Type: CR2032 (cod. **33212** optional)
- Dimension (LxWxH): 14x3,8x3 cm
- Weight: 70 g
- Norms: IEC60601-1, IEC60601-1-2 (EMC), ASTM E1965-98



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

