



ESF10 Line HEIGHT RODS



LITTLE ACCESSORIES



ESF10 Line HEIGHT RODS

Cod. ESF10-225
 < 9=, < HFC8

Measuring range	85 - 210 cm
Division	1 mm
Platform	280x320x45H mm



Cod. ESF10-222/F
 < 9=, < HFC8 fK 5@@HMD9L

Measuring range	92 - 203 cm
Division	1 mm



Cod. ESF10-208
 T9LESCOPI7 < 9=, < HFC8

Measuring range	0 - 200 cm
Division	1 mm



Cod. ESF10-209/F
 656M< 9=, < HFC8

Measuring range	0 - 100 cm
Division	1 mm



Cod. Z90012
 H5D9'A 95GI F 9

Measuring range	0 - 200 cm
Division	1 mm



Cod. ESF10-210/F
 656M< 9=, < HFC8

Measuring range	0 - 100 cm
Division	5 mm



ESSE3 srl, Via Garibaldi 30
 14022 Castelnuovo D.B. (AT)
 Tel +39 011 99 27 706
 Fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web: www.esse3-medical.com

