



ESB LINE – ACCESSORIES



ESB-688



ESB-689



ESB-689A



ESB-685



ESB-686



ESB-684



ESB-683



ESB-004



ESB-005



ESB-009



ESB-015

ESB-688 **Overbid Table**

Specification: 750x400x850/1100mm

Material: steel powder coated, wooden.

ESB-689 **Gas-spring over bed Table**

Specification: 780x390x800/1300mm

Material: Steel powder coated, wooden

ESB-689A **over bed Table**

Specification: 600x400x95mm

Material: Stainless steel

ESB-685 **Working stool**

ESB-686 **stool**

Specification: Height adjustable 480-650mm

Material: Stainless steel

ESB-684 **Footstool**

Specification: 520x430x500mm

Material: Stainless steel

ESB-683 **Footstool**

Specification: 520x430x500mm

Material: Stainless steel

ESB-004 **X-Ray film viewing screen**

--- single

--- double

--- Three

ESB-005 **Visual Acuity Light E**

ESB-009 **Gloves shelf**

Specification: 560x360x1420mm

Material: Stainless steel

ESB-015 **Leg support**

Specification: 1125x610x890mm

Material: Stainless steel