



MEDICAL ACCESSORIES

ESH-4101

Saline Stand

- Stainless Steel Tubular Pipe
- Stainless Steel Adjustable Rod with Two Hooks
- Five Legs MS base epoxy powder coated on 5cm Castors
- Option: **MF4101A**
Tubular Epoxy powder coated MS pipe & chrome plated Adjustable rod with two hooks



ESH-4102

Saline Stand

- Stainless Steel Tubular Pipe
- Stainless Steel Adjustable Rod with Two Hooks
- Five Legs MS base epoxy powder coated on 5cm Castors
- Option: **MF4102A**
Tubular Epoxy powder coated MS pipe & chrome plated adjustable rod with two hooks



ESH-4200

Bed Side Screen

(3 Panels)

- Size: 165 H x 240 L cms
- Tubular Frame work mounted on 5 cms Castors with Curtains
- Framework of Mild Steel, epoxy powder coated
- On Castors
- Supplied with 1 set of curtains



ESH-3201

Bed Side Screen

(4 Panels)

- Size : 165 H x 240 L cms
- Tubular Framework mounted on 5cms castors with Curtains
- Framework of Mild Steel, epoxy powder coated
- On Castors
- Supplied with 1 set of curtains





ESH-4300
Kick Bucket

- Frame work of Stainless Steel with Stainless Steel Bowl
- On three castors
- Option: **MF4300A** Framework of Mild Steel with epoxy powder coating



ESH-4301
Kick Bucket

- Framework of Stainless Steel with Stainless Steel Bucket
- On four castors
- Option: **MF4301A** Framework of MS with epoxy powder coating



ESH-4400
Foot Step Double

- Frame work of Stainless Steel
- Rubber Stumps



ESH-4401
Foot Step Double

- Frame work of Mild Steel, Epoxy Powder Coated
- Rubber Stumps
- Anti-Slip Platform



ESH-4402
Foot Step Single

- Size : 18" x 9" x 9"
- M.S. Tubular Legs Fitted with PVC Stump
- Pre Treated and Epoxy Powder Coated
- Anti-Slip Platform



ESH-4403
Foot Step Double

- Frame work of Mild Steel, Epoxy Powder Coated
- Rubber Stumps
- Anti-Slip Platform



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com

