



ESN Line

LED x-ray Film Viewers



MEDICAL FURNITURE



ESN Line

LED x-ray Film Viewers



Cod. **ESN-03643**



Cod. **ESN-07143**



Cod. **ESN-10643**



Cod. **ESN-14143**

Features:

- Devices in conformity with 2006/95 CEE (low voltage), and 2004/108 CEE (electromagnetic compatibility).
- Structure in light blue aluminium.
- Diffuser in white opal methacrylate, 3mm thickness.
- X-ray film support with stainless steel clips inox applied within the structure.
- Brightness Adjustment from 600 to 6000 lux.
- Display digit from 0 to 9 indicating the adjustment degree.
- Each click increases of 600lux.
- 7000°K color temperature.
- Storage of the adjustment degree when the unit is switched off.
- Automatic or manual power on/off from frontal controls for each sector, after turned on the main switch.
- Storage of the switching on choice with display indication.
- IP 5X protection.
- Isolation Class 1° - type B.

Electrical Components:

- Socket with fuse 10A – 250V
- Black bipolar switch D20 – 0-1 – 10A – 250V
- 2m black power cable – 3x0,75mm² – 10A – 250V
- LED bars SMD50x50
- Power supply AC170 - 250V – 12V – 30W – IP67
- Regulator board 30+30W
- Glass Fuses 5x20

Technical features:

Cod.		ESN-03643	ESN-07143	ESN-10643	ESN-14143
Light sectors		1	2	3	4
Light section size	cm	36 x 43	71 x 43	106 x 43	141 x 43
Total size	cm	47,3 x 54,2 x 4,8	82,3 x 54,2 x 4,8	117,3 x 54,2 x 4,8	152,3 x 54,2 x 4,8
Power supply	V/Hz	230/50	230/50	230/50	230/50
Power Consumption Min/max	W	16/67	21/134	26/192	30/255
Absorbed current	A	0,27	0,47	0,75	1,1
Fuses	Num/A	2/0,8	2/1	2/1,6	2/2
Led bars	Num/W	12/5,5	24/5,5	36/5,5	48/5,5
Power Adapters	Num/W	2/30	4/30	6/30	8/30



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

