



## Cod. ESF2-362/52


Epoxy painted structure  
 Tubular frame 50x25x1,5 mm with two adjustable feet Ø 50.  
 Plan of use with high density padding, Skai covered, washable with common detergents.  
 Two reclining sections with mechanical self-locking clutch system.  
 Backrest with an inclination of 55°.  
 The third section is foldable.  
 Paper roll support placed under the backrest.  
 Face hole (OPTIONAL).  
 Standard or anatomic knee crutches with sliding support and rods (OPTIONAL).



### DIMENSIONS

Length	185 cm
Width	60 cm
Height	85 cm
Capacity	Kg 130
Weight	36 Kg





## Cod. ESF2-362/52

Structure en tube d'acier émaillé.  
Châssis tubulaire de 50x25x1,5 mm avec deux pieds réglables de Ø 50.  
Plan d'utilisation avec un rembourrage à haute densité, revêtement en Skaï, lavable aussi avec les détergents plus communs.  
Deux sections inclinables avec système d'embrayage autobloquant mécanique.  
Dossier avec une inclinaison de 55°.  
La troisième section est pliable.  
Le support du rouleau de papier est placé sous le dossier.  
Trou pour le visage (EN OPTION).  
Appui-jambes standard or anatomiques avec support coulissant et tiges (EN OPTION).

### DIMENSIONS

Longueur	185 cm
Largeur	60 cm
Hauteur	85 cm
Capacité	Kg 130
Poids	36 Kg

