



EXAMINATION TABLES

ESH-3604

Gynaec Examination Table (Three Sections)

- Overall Size : 180L x 60W x 46H cms
- Framework made of Stainless Steel
- Back rest section adjustable on Ratchet / Lever
- Legs fitted with PVC Stumps
- 5 cms foam cushioned top covered with rexine
- Adjustable knee crutches
- Supplied with a pair of shoulder supports



ESH-3605

Gynaec Examination Table (Three Sections)

- Overall Size : 180L x 60W x 46H cms
- Framework made of Stainless Steel
- Back rest section adjustable on Ratchet / Lever
- Legs fitted with PVC Stumps
- 5 cms foam cushioned top covered with rexine with perineal cut
- Adjustable lithotomy poles with straps
- Finish : Epoxy Powder Coated



ESH-3506

Examination Table in 2 Parts (2 Section Top)

- Approx. Dimensions: 180L x 60W x 80H cms
- Adjustable backrest
- Two section cushioned top
- Legs fitted on PVC Stands
- Stainless Steel Framework



ESH-3507

Examination Table in 2 Parts (2 Section Top)

- Approx. Size : 180L x 60W x 80H cms
- Adjustable backrest
- Two section foam padded, upholstered top
- Pre treated & Epoxy powder coated
- Mild Steel Framework



ESH-3612

Obstetric Delivery Bed in 2 Parts

(2 Section Top)

- Approx. Dimensions: 72"L x 30"W x 30"H
- Tubular frame work made in two section
- Body section with legs mounted on PVC stumps
- Leg section on swivel castors with brakes
- Two section cushioned top
- Supplied with a pair of knee crutches/Leg Stirrups
- Mild Steel Framework
- Pre treated & Epoxy powder coated



ESH-3608

Gynaec Examination Table

(Two Section)

- Overall Size : 180L x 60W x 80H cms
- Framework made of rectangular / Square M.S. Tube
- Back rest section adjustable on Ratchet / Lever
- Legs fitted with PVC Stumps
- 5 cms foam cushioned top covered with rexine with perineal cut
- Adjustable lithotomy poles with straps
- Finish : Epoxy Powder Coated



ESH-3610

Examination couch with cabinet & drawers

(2 Section Top)

- Approx. Dimensions: 187L x 51W x 81H cms
- Adjustable backrest by gas spring system
- Two Section Foam Padded, Upholstered top
- Three Drawers, Three Cabinets
- Tray for B.P. Apparatus
- Pre treated & Epoxy powder coated



ESH-3611

Delivery Bed

(Fixed Height)

- Approx. Size –Body section -1000L x 800W x 750 H, Leg Section -900Lmm
- Bed, Labour and Delivery made in 2 section, body section with legs mounted on rubber tips
- Leg section swivel castors with brakes
- Leg section can be completely recessed under the body section supports for legs holders with solid mechanism to fix leg holders at required height & width
- Filled mattress – pad and cover mattress pad in 2 section
- Supplied with 2 knee crutches adjustable in height & width knockdown construction
- Epoxy powder coated finish



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com

