



ES -COMMUNE CHAIRS



ES609

Foldable chromed steel frame, soft commode seat, fixed armrest and footrest, solid castor and rear wheel



ES654LGC

Foldable aluminium frame, reclining high back, soft commode seat, detachable armrest, elevating footrest, solid castor, pneumatic rear wheel, solid commode seat



ES683L

Aluminium chair frame, soft commode seat, FLIP-UP desk armrest, detachable footrest, solid castor, pneumatic rear wheel



ES609L

Foldable aluminium chair frame, soft commode seat, fixed armrest and footrest, solid castor, pneumatic rear wheel

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
ES609	64	44	58	20	49	87	100	94x29x89	19	21,2	120	244
ES654LGC	64	46	61	20	52	122	100	94x41x86	17	19	90	184
ES683L	61	46	57	15	49	93	100	79x34x77	14,2	16,4	115	240
ES871	61	41	57	15	52	90	100	94x33x90	17	19	120	244



ES699L

Aluminium chair frame, fixed armrest foldable footrest, solid rear castors with lock, plastic commode seat



ES696

Foldable chromed steel frame, fixed armrest, solid rear castors with lock, plastic commode seat, adjustable seat height



ES692

Powder coating steel frame, FLIP-DOWN armrest, detachable footrest, solid rear castors with lock, plastic commode seat



ES691

Powder coating steel frame, detachable armrest and footrest, solid rear castors with lock, plastic commode seat

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
ES699L	54	45	7,5	7,5	49	91,5	100	56x56x72	14,5	16	120	240
ES809B	55	45	7,5	7,5	49-59	91-101	100	52x23x76	7,3	8,4	200	400
ES870ABJ	48	45	12,5	12,5	53	96	100	56x54x58	15	17	160	320
ES868	48	44	12,5	12,5	52	94	100	54x53x53	17	19	160	320



ES693

Powder coating steel frame , fixed armrest , foldable footrest , solid rear castors with lock , plastic commode seat



ES695S

Stainless steel chair frame , plastic commode seat FLIP-DOWN armrest , foldable footrest , detachable back-rest , rear castors with lock ,



ES893

Chromed steel frame , plastic commode seat with detachable cover , fixed armrest



ES895L

Aluminium chair frame, plastic commode seat with detachable cover , fixed armrest , adjustable height

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
ES693	55	44	12,5	12,5	48	84	100	71x29x67	12	13	180	378
ES695S	48	44	12,5	12,5	55	102	100	97x31x85	16,1	18,3	88	192
ES893	48	44			47	84	100	50x19x75	8	8,5		
ES895L	48	48			47-51	70	100	60x38x52	5	5,5		





ES - SHOWER CHAIRS



ES792



ES7920L



ES7982L



ES795



ES796L



ES7931



ES799L



ES790



ES791S



ES793S



ES794L



ES7962L



ES798L



ES797L



ES – WALKING AIDS



ES9632L



ES963L



ES915L



ES919L



ES913L



ES976L



ES9125L

ES912L



ES910L

ES920L

ES930L

ES929L

ES939L

ES923L

ES937L

ES9300L

ES933L

ES948L



ES945

ES932

ES926

ES942S

ES922

ES946S



ESSE 3 Via Garibaldi 30
 14022Castelnuovo D.B. (AT)
 tel +39 011 99 27 706
 fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web : www.esse3-medical.com

