



ES – SPORT WHEEL CHAIRS



ES730LQ-36



ES730LQ-36



ES873



ES871

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
ES730LQ-36	68	36	61	10	48	77	100	76x62x54	9,6	11,6	108	230
ES720LQ-36	66	36	59	10	48	72	100	87x36x73	10,8	13	120	240
ES21LQ-36	60	36	59	10	48	72	100	87x36x73	12,1	14,3	120	240
ES730LQ-36	60	36	59	10	44	72	100	87x36x73	12,1	14,3	120	240



ES731LQ-36



ES723LQ-36

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
ES731LQ-36	68	36	61	10	52	84	100	70x61x54	10	12	108	230
ES732LQ-36	60	36	61	10	48	78	100	87x36x79	11,7	13,9	120	240



Shooting



ES744LQ-43



Fencing

ES766LQ-36



ES740LO-36



ES756LQ-36



ES755LQ-36



ES785LQ-36



ES778LQ-36



Basketball Guard



Basketball Center

ES779LQ-36

ES777LQ-36


Rugby Offensive



ES770LQ-32


Rugby Offensive



ES771LQ-32



ES786LQ-36



ES710L-30



ES – POWER WHEELCHAIRS



ES106L

Powder coated aluminium frame , detachable armrest , pneumatic castor , pneumatic drive wheel, detachable footrest , lithium battery , foldable chair frame . Max load 110 kg



ES107LA

Powder coated aluminium frame , detachable inclined armrest , detachable footrest , pneumatic castor , pneumatic drive wheel, lithium battery , foldable chair frame . Max load 110 kg



ES108LA

Powder coated aluminium frame , detachable inclined armrest , detachable footrest , pneumatic castor , pneumatic drive wheel, lithium battery , foldable chair frame . Max load 110 kg



ES101A

Powder coated steel frame , detachable armrest , detachable footrest , rubber foam castor , inflat-free hollow drive wheel with hand rim . Drop back handle , foldable chair frame . Max load 110 kg

MODEL N°							SPEED km/h	BATTERY	W	PACKAGE CM	NW KG	GW KG	CTNS 20" FCL	CNTS 40" FCL
ES106L	60	41	61	17	49	76	7	20AH	100x02	106x36x78	19,7	21,9	64	198
ES907BC	63	43	61	17	50	92	7	20AH	100x2	88x36x71	17,2	19,4	96	234
ES909	63	43	61	17	50	92	7	20AH	100x2	88x36x71	17,9	20,1	130	310
ES908AQ	60	46	50	20	50	92	7	30AH	200x2	89x48x82	16,9	19,1	44	88



Autonomy 20 km

Max slope 12°

ES121

Powder coated steel frame, FLIP-UP desk armrest, detachable footrest, rubber foam castor and drive rear wheel, foldable backrest. Max load 130 kg

ES111A

Powder coated steel frame, drop back handle, FLIP-UP armrest, detachable footrest, rubber foam castor and drive rear wheel, foldable chair frame. Max load 120 kg



Autonomy 25 km



Max slope 12°

ES123GC-43

Powder coated steel frame, adjustable armrest PAD, detachable sideboard and footrest, angle adjustable backrest, rubber foam castor and drive rear wheel, with safety lamp. Max load 130 kg

ES123-43

Powder coated steel frame, adjustable armrest PAD, detachable sideboard, detachable and elevating footrest, adjustable footplates, angle adjustable backrest, pneumatic castor and drive rear wheel, adjustable headrest, with safety lamp. Max load 120 kg

MODEL N°							SPEED km/h	BATTERY	W	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
ES121	63	46	31	20	53	94	6	28AH	200x2	83x63x74	49	579	46	96
ES111A	60	41	31	20	50	91	6	28AH	200x2	85x44x81	58,0	66,6	50	104
ES123GC-43	67	43	35	22	58	100	8	38AH	200x2	90x74x100	75	85	28	90
ES123-43	65	43	31	20	52	99	6	28AH	200x2	82x69x79	62	70	81	216



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.YggY' la YXJWJ.com

