



# ES -WHEEL CHAIRS



**ES802-35**

Children type, chromed steel , fixed armrest , fixed footrest , solid castors, solid rear wheel



**ES875**

Powder coating steel frame , fixed armrest , fixed footrest , solid castors, solid rear wheel



**ES873**

Chromed steel frame , Powder coating cross bar , fixed armrest , fixed footrest , solid castors, solid rear wheel



**ES871**

Chromed steel frame , fixed armrest , fixed footrest , solid castors, solid rear wheel

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
<b>ES802-35</b>	55	35	55	15	47	87	75	88x23x89	16,2	18,4	164	410
<b>ES875</b>	60	41	59	20	50	87	100	94x25x89	17,3	19,5	132	336
<b>ES873</b>	60	41	59	20	50	87	100	94x25x89	17,3	19,5	144	372
<b>ES871</b>	60	41	59	20	50	87	100	94x25x89	17,3	19,5	144	372



**ES869X**

Powder coating steel frame, double cross bar, fixed inclined armrest, fixed footrest, solid castors, solid rear wheel



**ES873**

Chromed steel frame, fixed armrest, fixed footrest, solid castors, solid MAG rear wheel



**ES870ABJ**

Chromed steel frame, double cross bar, drop back handle with brake, fixed armrest, fixed footrest, solid castors, pneumatic MAG rear wheel



**ES868**

Powder coating steel frame, double cross bar, fixed armrest, fixed footrest, solid castors, solid rear wheel

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
<b>ES869X</b>	61	41	59	20	49	87	100	94x23x89	17	19	144	372
<b>ES809B</b>	60	41	61	20	50	88	100	94x25x89	18,7	20,9	140	336
<b>ES870ABJ</b>	60	41	40	15	45	89	100	82x28x76	11	13	168	387
<b>ES868</b>	59	46	59	20	47	91	100	94x27x93	16	18	108	272



**ES874**

Chromed steel frame , double cross bar , fixed armrest , fixed footrest , solid castors, solid rear wheel



**ES876**

Chromed steel frame , double cross bar , fixed armrest , fixed footrest , pneumatic front wheel , solid rear castors



**ES874B**

Chromed steel frame , double cross bar , fixed armrest , fixed footrest , solid castors, solid MAG wheel



**ES874A**

Chromed steel frame , double cross bar , drop back handle , fixed armrest , fixed footrest , solid castors, solid rear wheel

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
<b>ES874</b>	60	41	59	20	47	91	100	94x28x92	17,4	19,6	108	272
<b>ES876</b>	62	39	20	61	45	84	100	97x31x85	16,1	18,3	88	192
<b>ES874B</b>	61	41	61	20	47	91	100	94x28x92	18,7	20,9	108	272
<b>ES874A</b>	60	41	59	20	47	91	100	94x28x78	17,4	19,6	108	336



**ES909B**

Powder coated steel frame , FLIP-UP armrest , detachable footrest , solid castors, solid MAG wheel



**ES907BC**

Chromed steel frame , detachable desk armrest, elevating footrest , solid castors, pneumatic MAG wheel



**ES909**

Powder coated steel frame , FLIP-UP desk armrest , detachable footrest , solid castors, solid rear wheel



**ES908AQ**

Chromed steel frame , FLIP-UP desk armrest , detachable footrest , solid castors, quick release pneumatic rear wheel, drop back handle

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
<b>ES909B</b>	57	41	51	20	50	93	100	82x28x95	19,7	21,9	130	310
<b>ES907BC</b>	58	41	40	20	50	91	100	78x28x91	17,2	19,4	180	366
<b>ES909</b>	56	41	59	20	50	90	100	83x28x95	17,9	20,1	130	310
<b>ES908AQ</b>	58	41	61	15	50	93	100	80x28x78	16,9	19,1	168	402



**ES951B-56**

Powder coated steel frame , 22' seat width , FLIP-UP armrest , detachable footrest , rubber foam castors, pneumatic MAG wheel



**ES975-51**

Chromed steel frame , 20' seat width , detachable armrest , detachable footrest , solid castors, pneumatic rear wheel



**ES210ABE-61**

Powder coated steel frame , detachable desk armrest , detachable footrest , 24' seat width , foldable push bar , solid castors, pneumatic MAG wheel



**ES209AE-61**

Powder coated steel frame , detachable desk armrest , detachable footrest , 24' seat width , foldable push bar , solid castors, pneumatic rear wheel

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
<b>ES951B-56</b>	71	56	60	20	54	94	125	82x35x95	22,6	24,8	112	245
<b>ES975-51</b>	70	51	60	20	50	92	125	80x28x93	20	22	130	310
<b>ES210ABE-61</b>	78	61	30	20	50	96	125	73x32x83	19,3	21,5	144	324
<b>ES209AE-61</b>	80	61	61	20	50	96	125	86x37x83	22,2	24,4	100	240



**ES958LBHP**

Aluminium chair frame , reclining high back , adjustable headrest , detachable armrest , elevating footrest , solid castors, pneumatic MAG wheel , adjustable seat angle, safety belt

**ES203BJP**

Powder coated steel frame , reclining high back , adjustable and detachable headrest , detachable pillow , adjustable seat angle , height adjusting armrest , elevating footrest , solid castors, pneumatic MAG wheel , united brake



**ES985LBJ**

Aluminium chair frame , reclining high back, adjustable headrest , detachable armrest , fixed footrest , adjustable seat angle, solid castors, pneumatic MAG wheel , united brake



MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
<b>ES958LBHP</b>	52	32	40	15	50	118	100	85x39x95	20,8	23	90	223
<b>ES203BJP</b>	63	46	61	20	53	118	150	82x67x84	35,3	37,7	70	148
<b>ES985LKBJ</b>	52	38	40	15	50	118	75	85x39x89	18,3	20,5	90	223



**ES980LA-35**

Children type aluminium chair frame , FLIP-UP armrest , fixed footrest , solid castors, pneumatic rear wheel , drop back handle



**ES980LA**

Aluminium chair frame , FLIP-UP armrest , fixed footrest , solid castors, pneumatic rear wheel , drop back handle



**ES801LB**

Mini foldable portable aluminium chair frame , FLIP-UP armrest , solid castors, pedal locking rear MAG wheel , with seat belt , leg straps



**ES803LB**

Mini foldable portable aluminium chair frame , FLIP-UP armrest , solid castors, pedal locking rear MAG wheel , with seat belt , leg straps

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
<b>ES980LA-35</b>	51	35	57	15	48	89	75	90x33x73	13,6	16	148	315
<b>ES980LA</b>	57	41	61	15	50	91	100	94x33x73	15	17,2	126	284
<b>ES801LB</b>	38	32	12,5	10	45	82	100	58x39x26	6,2	8,2	496	1200
<b>ES803LB</b>	38	32	15	12,5	45	84	100	71x39x27	6	7,5	357	148



**ES869LX**

Aluminium light weight chair frame , fixed inclined armrest, solid castors, pneumatic rear wheel , fixed footrest



**ES8864L**

Aluminium light weight chair frame , fixed armrest, solid castors, pneumatic rear wheel ,



**ES871LB**

Aluminium light weight chair frame , fixed armrest, fixed footrest , solid castors, pneumatic MAG wheel



**ES871LBJ**

Aluminium light weight chair frame , fixed armrest, fixed footrest , solid castors, pneumatic MAG wheel with united brake

MODEL N°							MAX LOAD kg	PACKAGE CM	NW KG	GW KG	CTNS 20' FCL	CNTS 40' FCL
<b>ES869LX</b>	61	41	59	20	49	87	100	94x23x89	13	15	144	372
<b>ES801LB</b>	61	41	60	20	45	89	100	94x29x92	12,7	14,9	108	272
<b>ES871LB</b>	58	41	40	15	47	88	100	80x28x90	10	12	140	340
<b>ES871LBJ</b>	58	41	40	15	47	89	100	80x28x90	12	15	140	340