



NEBULIZER

FOR INTENSIVE USE

Cod. F-205



COMPRESSOR NEBULIZER FOR AEROSOLTHERAPY FOR INTENSIVE USE

It is equipped with piston type compressor. Powerful, silent and reliable.

ABS case with large built-in accessories compartment.
Supplied with accessory kit:
nebulizer, adult and child mask,
210 cm tubing.

Cod.	F-205	F-205 PRO
Dimensions	300 x 190 x 160 (H) mm	300 x 190 x 160 (H) mm
Weight	2,1 kg	2,3 kg
Power supply	230V - 50/60 Hz	230V-50Hz • 230V-50/60Hz
Power consumption	185 VA	195 VA
Nebulization	2 cc - 7 minute	3 cc - 7 minute
Operating cycle	CONTINUOUS	45 min
Working pressure	2,5 bar	3,8 bar(50Hz) / 3 bar(50-60Hz)

Optional Accessories



- 07 NEBULIZER BULB
- 07 MASK ADULT
- 07 MASK CHILD
- 07 TUBE 210 cm
- 07 FILTER AEROSOL



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com





NÉBULISEUR À USAGE INTENSIF

Cod. F-205



NÉBULISEUR POUR AÉROSOLTHÉRAPIE À USAGE INTENSIF

Equipé de compresseur à piston,
puissant, silencieux et de longue
durée.

Coffret ABS avec compartiment
pour accessoires.

Accessoires inclus: ampoule,
masque adulte et enfant, tuyau de
210 cm

Cod.	F-205	F-205 PRO
Dimensions	300 x 190 x 160 (H) mm	300 x 190 x 160 (H) mm
Poids	2,1 kg	2,3 kg
Alimentation	230V - 50/60 Hz	230V-50Hz • 230V-50/60Hz
Puissance absorbée	185 VA	195 VA
Nébulisation	2 cc - 7 minute	3 cc - 7 minute
Cycle d'opération	CONTINUOUS	45 min
Pression maximale	2,5 bar	3,8 bar(50Hz) / 3 bar(50-60Hz)

Accessoires en Option



- 07 NEBULIZER BULB
- 07 MASK ADULT
- 07 MASK CHILD
- 07 TUBE 210 cm
- 07 FILTER AEROSOL



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

