



9G!5I !' \$@
6 YbW 'hcd' ghYUa '
ghYf]]nYf



STERILIZATION



9G!5I !' \$@

''''6YbW `hcd `ghYUa `ghYf`]nYf''''



- LCD display, 8 programs
- Automatic feeding water and draining water system
- Full-Automatic operation and door lock system
- Printer and USB port
- Automatic and safety door lock system
- According to Europe B standard, with 3 times pre-vacuum
- Failure reason displayed on LCD
- New type of separate steam generator
- Multi-language
- Detailed print with monitoring temperature and pressure record in 6 seconds.

Technical Features

Chamber Dimensions	250×600 mm (Φ×D)
Chamber Capacity	30L
Overall Size	48×80×42 cm (L×W×H)
Heater	2000W
Net Weight	55 kg
Gross Weight	65 kg
Packing Size	94×57×50 cm
Minimum Desk Size	60×50 cm
Power Supply	220±10%,50HZ
Sterilize Temperature	121/134
Sterilization Time	4-18mins
Dry Procedure	Vacuum drying
Sterilizer Construction	Stainless steel 304
Clean Water tank	Internal
Stainless steel instrument container	Strays
Safety	<ul style="list-style-type: none"> • Over heat control: When inner temperature is overheated, power will be shut down automatically with an alarm system • Over pressure control: alarm in case of over pressure. • Safety door lock system • Low water level indicator: alarm in case of low water



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

