



ESSELYA IFE

Electrophoresis Automated System



LABORATORY



Versatile

Allows to perform the most common techniques: Serum proteins, Lipoproteins, Urinary Proteins, Immunofixation.

Complete Automation

ESSELYA FE automatically arranges the loading of all the reagents, to standardize the delicate electrophoretic process.



Serum Proteins



Haemoglobins



Lipoproteins



Integrated webcam
for reporting of results



Automated Ag/Ab
reaction contact

Productivity

Serum-Proteins : 40 Samples

Lipoproteins: 40 Samples

Haemoglobins: 40 Samples

Immunofixation: 5 Samples



Software features

User interface: monitor, keyboard and mouse

Minima recognition: automatic

Managing techniques: adjustable techniques in times,
migration voltage/current

Operating system:



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

