



ELECTRIC OPERATING TABLE

FOR GYNECOLOGY AND OBSTETRICS

Cod. **ESK-A99/5**



TECHNICAL FEATURES

Electric operating table gynecology and obstetrics.

The table top elevation, trendelenburg/anti trendelenburg and backrest are all controlled by both micro touch remote control and foot switch. Two electric motors that enable flexibility, smooth movements, low noise, and high reliability. All the base and column are covered in 304 medical grade stainless steel.

Mattress with different available colors, easy to clean and sterilize. Third section extensible.

Length: 2000 mm

Width: 600 mm

Elevation (up/down): 920mm/650mm

Trendelenburg/anti-trendelenburg: 8° / 20°

Backrest (up/down): 45°/7°

Accessories :

- 1 Mattress
- 2 Armrests
- 2 Adjustable knee crutches
- 1 Leg Support
- 1 Removable stainless steel tray
- 2 Handgrips
- 4 Wheels, 2 with brakes
- 1 Footswitch
- 2 Pedals
- 1 Remote Control



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com





TABLE D'OPÉRATION ÉLECTRIQUE

POUR LA GYNECOLOGIE ET L'OBSTETRIE

Cod. ESK-A99/5



CARACTERISTIQUES TECHNIQUES

Table d'opération électrique pour la gynécologie et obstétrique.

L'élévation de la table, trendelenburg / Anti trendelenburg et le dossier sont tous contrôlés à la fois par micro contact télécommande et l'interrupteur au pied. Deux moteurs électriques qui permettent la flexibilité, des mouvements lisses, faible bruit, et une grande fiabilité. La base et la colonne sont recouvertes d'acier inoxydable 304 de qualité médicale.

Matelas disponible en différentes couleurs. Troisième section extensible.

Longueur: 2000 mm

Largeur: 600 mm

Élévation (en haut/bas): 920mm/650mm

Trendelenburg/anti-trendelenburg: 8° / 20°

Élévation du dossier (en haut/bas): 45°/7°

Accessoires :

- 1 Matelas
- 2 Appuie-bras
- 2 Appuie-jambes anatomiques réglables.
- 1 Support pour les jambes
- 1 Plateau amovible en Acier inoxydable
- 2 Poignées
- 4 Roues dont 2 avec frein
- 1 Commutateur au pied
- 2 Pédales
- 1 Télécommande



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

