



ES-OM

Oxygen Monitor



INTENSIVE CARES



ES-CA

Oxygen Monitor

The ES-OM is a device applicable for measuring oxygen concentration and gas level in various places, even in infant incubator. Compact size & easy to use.



Technical Specifications:

Cod. ES-OM



Measure Range	0-100% O ₂ in air
Response Time	T ₉₀ 13 seconds in 25
Repeatability	1%
Linear Error	3%
Operation Temperature	0-45°C
Display	STN Blue LCD with CCFL back light
Display area	120 x 90 mm
Resolution	320 x 240 pixels
Writing speed	25mm/sec
Alarms	Visual (Bell) and Audio (Buzzer beep)
Key board	Flat touch membrane operation keys
Battery life	more than 24 hours
Battery	3AA alkaline battery or NH rechargeable battery
Dimensions	235 x 166 x 140mm



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

