



FINGERTIP OXIMETER



Cod. 34282



**With
LCD Backlit
Display**

In order to monitor pulse oxygen saturation and pulse rate through patient's finger (adult/child over 5 years). Small, light, easy to use it is ideal in hospital, medical center, home, before and after sport.

- Backlit display
- Lower voltage alarm display with one-key operation
- High accuracy and repeatability
- Automatic power-on and power-off in 8 seconds (if finger not in)
- SpO2 and pulse monitoring, and histogram of pulse display
- Audible and visible alarms
- P.I. perfusion index
- Pediatric version cod. 34265

TECHNICAL SPECIFICATIONS	
Oxygen saturation range	50% to 99%
Pulse rate range	30 bpm to 240 bpm
Backlit LCD Display	oxygen saturation, pulse rate, low battery indicator
SpO2 accuracy	±2% between 75% to 99% ±3% between 50% to 75%
Pulse accuracy	± 2 bpm
Power	2 AAA 1.5 V alkaline batteries (for 100 hours of continuous operation)
Dimensions	6,6 x 3,6 x 3,3 cm
Weight	60g



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com





OXYMÈTRE DE POULS

Cod. **34282**



Avec
Écran LCD
Rétroéclairé

Pour mesurer la saturation en oxygène et la fréquence cardiaque sur le doigt du patient adulte. Compact, léger, facile à utiliser, il est idéal dans les hôpitaux, les centres médicaux, à la maison, avant et après le sport.

- Mesure de la SpO₂ et la fréquence cardiaque, et affichage de l'histogramme du pouls
- Histogramme de la fréquence cardiaque
- Calcul de l'index de perfusion
- Alarmes vidéo et audio
- Hautes précision et reproductibilité
- Extinction automatique en 8 secondes (après retrait du doigt)
- Mesure de la SpO₂ et la fréquence cardiaque, et affichage de l'histogramme du pouls
- Version pédiatrique avec code **34265**

SPECIFICATIONS TECHNIQUES

Plage de saturation en oxygène	50% à 99%
Plage de fréquence cardiaque	30 bpm à 240 bpm
Affichage LCD	saturation en oxygène, fréquence cardiaque, indicateur de batterie faible
Précision SpO ₂	±2% entre 75% et 99% ±3% entre 50% et 75%
Précision fréq. cardiaque	± 2 bpm
Alimentation	2 piles alcalines AAA 1,5 V (pour 100 heures de fonctionnement continu)
Dimensions	6,6 x 3,6 x 3,3 cm
Poids	60g



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

