



# ES-BS5

Handheld digital  
medical scope



LITTLE ACCESSORIES



# ...9G!6G)

## Handheld digital medical scope

### Product Features

Specially designed for the department of ENT. It can be used to examine eardrum, nasal cavity and throat with different camera probes.

- Avoid high cost of the traditional ENT examination equipment, such as light source, CCD, Monitor etc. All the functions are in one handheld device.
- Light weight and portable, suitable for outpatient department, emergency room and medical ward.
- High resolution image. The SD card can store multiple frozen images and videos, and can be connected to computer to print patient reports.
- 3 inches high definition LCD screen, watching and taking photos at the same time.
- It can be connected to larger size monitor, synchronized to the external bigger screen.

### Product Specifications

- Charge: 11 OV-220V AC 50/60HZ Output DC 5V 500MA
- LCD: 3 inches special LCD screen, minimum color difference for different viewing angle.

### Applicable Camera Probes

- The device is applicable for viewable examination for ear, nose and throat with different camera.



Otoscope (standard)



Rhinoscope (optional)



Laryngoscope (optional)



Cod. **ES-BS5**



Normal ear drum



Normal ear drum



Ear drum with otitis media



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

