



GLUCOSE MONITOR

ON CALL PLUS



Advanced technology glucose monitors. Available as unit only or complete kit.

The kit includes 10 strips, code chip, control solution, 10 sterile lancets, lancet device, clear cap (for testing on forearm or palm), carrying case, user manual (English & French), logbook.

Main features:

- accurate result in just 10 seconds
- only 1 μ l sample size
- automatic calibration with code chip
- plasma calibration results
- advanced data management with 7, 14, 30 days average
- affordable for frequent monitoring
- flexibility to test on the forearm, palm or fingertip

Technical Specifications

- Measurement Range: 20 to 600 mg/dL or 1.1 to 33.3 mmol/L
- Sample Sites - Size: Finger tips, Forearm or palm - 1 μ L
- Test Time. Approximately 10 seconds
- Power Source: One 2032 or equivalent 3.0 V coin cell battery
- Battery Life: 12 months or approximately 1.000 tests
- Memory Plus: Up to 300 records with time and date



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com





CHOLESTEROL, TRIGLYCERIDES, AND GLUCOSE MONITOR **MULTICARE**



Unique diagnostic device for self-check of glucose, cholesterol and triglycerides. Multi-parameter technology allows to measure, with the same reflectometer, different parameters by simply changing type of strips: a chip in each pack of strips is specific for each analysis. Regular measuring of these three parameters is important to prevent cardiovascular pathologies if we consider that intolerance to glucose, hyper-cholesterolemia and hyper-triglyceridemia are three of the main cardiovascular risk factors, along with age, obesity, lifestyle, smoking and arterial hypertension.

Multicare is supplied in 3 language carton box (English, Italian, Spanish) with a nylon bag, lancet device, 10 lancets, user manual (English, Italian, Spanish as above). Available on request French and Portuguese manual.

Technical Specifications

- Measuring range:
 - glucose: 25-500mg/dl (1.3-28 mmol/l)
 - triglycerides 50-500 mg/dL (0.5-5.6 mmol/l)
 - cholesterol 130-400 mg/dL (3.3-10.2 mmol/l)
- LCD display
- Memory: 255 test results with date and time
- Size: 48 x 104 x 20.5 mm
- Weight: 65 g
- Batteries: 2 batteries



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

