




# ESF1-150/002

**BLOOD DONOR AND RELAX CHAIR**



**Cod. ESF1-150/002**

## BLOOD DONOR AND RELAX CHAIR

Made of washable (class 1 IM) and disinfectant treatable material.

Anatomical padding with different density.

Adjustable armrests height can rotate 180°.

Tilt in Trendelenburg and reverse Trendelenburg position:

for backrest between +60° and -5°

and for legs section between -80° and +5°

Maximum load capacity: 130 Kg.

Dimensions: 55 x 175 x 51(H) cm




ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)





# ESF1-150/002

**FAUTEUIL PRÉLÈVEMENT ET RELAX**



**Cod. ESF1-150/002**

## FAUTEUIL PRÉLÈVEMENT ET RELAX

Fabriqué en matériaux lavables, (classe 1 IM).

Fauteuil rembourré à densités différentes et anatomiques.

Les appui-bras sont réglables en hauteur et ils peuvent pivoter de 180°.

Inclinaison en position Trendelenburg / anti-Trendelenburg :

dossier réglable entre +60° et -5°,

repose pied réglable entre -80° et +5°. Capacité maximal : 130 kg

Capacité de charge maximale: 130 kg.

Dimensions: 55 x 175 x 51(H) cm



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

