



SB10 Projector

Specifications

Projection distance:	Adjustable up to 8 m
Magnification:	28x (in 5m refraction)
Projection field:	210X308mm (in 5m refraction)
Screen size:	380X420mm
Angle of inclination:	±15°
Number of test charts:	25
Masks:	Red/Green, polarized, single horizontal line, single vertical line, single character
Automatic standby	After 10 min., without operation
Bulb:	12V/30W Halogen
Size (Main unit):	193(H)×198(W)×263(L)mm
Weight (Main unit):	5.2Kg
Power Supply	AC220V/50Hz or AC110V/60Hz
Power consumption :	75VA (Max)
Infrared remote control:	39 Functions: 25 individual test options, 2 filter options, masks options, program test sequence, etc.

