



# AUTO-REFRACTOMETER

Cod. **ES-V0V815**



## Specifications:

<b>Measuring Range</b>	<b>Sphere: -20 ~ +20D (VD=12) 0,125D steps</b> <b>Cylinder: -8 ~ +8D 0,12D / 0,25D steps</b> <b>Axis: 0 ~ 180° 1°steps</b>
<b>Pupil Distance</b>	<b>45~85mm,1mm steps</b>
<b>Vertex Distance</b>	<b>0 mm, 12 mm, 13.75 mm</b>
<b>Min. Pupil Size</b>	<b>Ø 2,0 mm</b>
<b>Chart</b>	<b>Auto focus</b>
<b>Display</b>	<b>7"color LCD</b>
<b>Print</b>	<b>Thermal Printer</b>
<b>Power Supply</b>	<b>100 ~ 120v/200 ~ 240v,50/60Hz</b>
<b>Dimensions</b>	<b>484 (L) x 282 (W) x 440 (H) mm</b>
<b>Weight</b>	<b>~20kg (44 lb)</b>




**ESSE3** srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)





# AUTO-RÉFRACTOMÈTRE

Cod. **ES-VOV815**



## Spécifications:

Plage de mesure	Sphère: -20 ~ +20D (VD=12) 0,125D étapes Cylindre: -8 ~ +8D 0,12D / 0,25D étapes Axe: 0 ~ 180° 1° étape
Distance pupilles	45~85mm, 1mm étapes
Distance Vertex	0 mm, 12 mm, 13.75 mm
Taille minimum de la pupille	Ø 2,0 mm
Graphique	Mise au point automatique
Écran	7" color LCD
Imprimer	Imprimante thermique
Alimentation	100 ~ 120v/200 ~ 240v, 50/60Hz
Dimensions	484 (L) x 282 (P) x 440 (H) mm
Poids	~20kg (44 lb)



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

