



ODM2200



ULTRASONIC A/B SCAN



Cod. ODM2200

- 8 Color Codes
- High Resolution Images
- Operating Mode: B, B+A, B+B, A
- 8 Images Memory
- 6 Steps Zoom from 34-60mm
- Variable displaying Area
- Cineloop for 5.6sec
- Accurate Axial Measurement under B-mode (SMI-5)
- Manual Correction of Axial measurement under A-mode (SMI-5)
- Multi-groups of Calipers for Distance Measurement
- Various area measurement
- Measuring Mode: Contact and Immersion
- Save up to 50 groups of IOL calculation results
- 8 Groups' average with Standard Deviation (S.D.)
- Large Screen LCD Monitor
- B-Probe protection: Automatic stop after 10 minutes of continuous scanning
- Trackball for Easy Measurement and Operation
- Real-time Clock
- Optional Image Output: Video Printer or Workstation

PACHYMETER

Description :

Portable and sleek
High light large scale LCD display
User friendly
Auto gain control
Auto/manual mode
Built-in thermal printer
Data and Graph printout



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com






ODM2200



ULTRA SONOGRAFIE A/B SCAN



Cod. ODM2200

- 8 Codes Couleur
- Images Haute Résolution
- Mode de fonctionnement: B, B+A, B+B, A
- Mémoire 8 Images
- 6 Mesures du Zoom 34-60mm
- Zone affichant variable
- Cineloop pour 5.6 sec
- Mesure axial précise sous mode-B (SMI-5)
- Correction Manuel de Mesure Axial sous mode-A (SMI-5)
- Multi-groups de Calibres pour la Mesure de la Distance
- Mesure de surface Variable
- Mode de Mesure: Contact et Immersion
- Enregistrement de jusqu'à 50 groupes de résultats de calcul IOL
- 8 Groups' moyen avec Deviation Standard (D.S.)
- Moniteur LCD Grand écran
- Protection B-Probe: arrêt automatique après 10 minutes de balayage continu
- Trackball pour Mesure et fonctionnement faciles
- Horloge en temps réel
- Sortie de l'Image en option: Imprimante vidéo ou station de travail

PACHYMETER

Description :

Portable et élégant
Ecran LCD à échelle d'éclairage élevée
Contrôle de gain automatique
Mode automatique/manuel
Imprimante thermique incorporé
Impression graphique des données
Facile à utiliser



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

